

Product Datasheet

Human PDGFRL ELISA Kit (Ready to Use) (orb551959)

Description	Human Platelet Derived Growth Factor Receptor Like Protein (PDGFRL) Sandwich ELISA KitReady to Use: ...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q15198
Sensitivity	0.041ng/mL
Entrez	5157
Expiration Date	Please enquire.