

Product Datasheet

Human PARP4 ELISA Kit (Ready to Use) (orb551980)

Description	Human Poly ADP Ribose Polymerase 4 (PARP4) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9UKK3
Sensitivity	0.061ng/mL
Entrez	143
Expiration Date	Please enquire.