

Product Datasheet

Human PAOX ELISA Kit (Ready to Use) (orb551982)

Description	Human Polyamine Oxidase (PAOX) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for co...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	31.5pg/mL
Expiration Date	Please enquire.