

Product Datasheet

Human PCYOX1 ELISA Kit (Ready to Use) (orb552003)

Description	Human Prenylcysteine Oxidase 1 (PCYOX1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UHG3
Sensitivity	0.115ng/mL
Entrez	51449
Expiration Date	Please enquire.