

Product Datasheet

Human PIIICP ELISA Kit (Ready to Use) (orb552018)

Description	Human Procollagen III C-Terminal Propeptide (PIIICP) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	5.3pg/mL
Expiration Date	Please enquire.