

Product Datasheet

Human PDCD1LG2 ELISA Kit (Ready to Use) (orb552031)

Description	Human Programmed Cell Death Protein 1 Ligand 2 (PDCD1LG2) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9BQ51
Sensitivity	0.063ng/mL
Entrez	80380
Expiration Date	Please enquire.