

Product Datasheet

Human PDIA3 ELISA Kit (Ready to Use) (orb552096)

Description	Human Protein Disulfide Isomerase A3 (PDIA3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection r...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P30101
Sensitivity	0.055ng/mL
Entrez	2923
Expiration Date	Please enquire.