

Product Datasheet

Human PTPN6 ELISA Kit (Ready to Use) (orb552130)

Description	Human Protein Tyrosine Phosphatase, Non Receptor Type 6 (PTPN6) Sandwich ELISA KitReady to Use: Pre-...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P29350
Sensitivity	0.105ng/mL
Entrez	5777
Expiration Date	Please enquire.