

Product Datasheet

Human PKM2 ELISA Kit (Ready to Use) (orb552155)

Description	Human Pyruvate Kinase, Muscle (PKM2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P14618
Sensitivity	29pg/mL
Entrez	5315
Expiration Date	Please enquire.