

Product Datasheet

Human RAD23B ELISA Kit (Ready to Use) (orb552169)

Description	Human RAD23 Homolog B (RAD23B) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P54727
Sensitivity	0.052ng/mL
Entrez	5887
Expiration Date	Please enquire.