

Product Datasheet

Human RIPK3 ELISA Kit (orb552189)

Description	Human Receptor Interacting Serine Threonine Kinase 3 (RIPK3) Sandwich ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	11.7pg/mL
Expiration Date	Please enquire.