

Product Datasheet

Human RIPK3 ELISA Kit (Ready to Use) (orb552190)

Description	Human Receptor Interacting Serine Threonine Kinase 3 (RIPK3) Sandwich ELISA KitReady to Use: Pre-dil...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	11.7pg/mL
Expiration Date	Please enquire.