

Product Datasheet

Human RLN ELISA Kit (Ready to Use) (orb552202)

Description	Human Relaxin (RLN) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient an...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04808
Sensitivity	9.66pg/mL
Entrez	6013
Expiration Date	Please enquire.