

Product Datasheet

Human RXFP3 ELISA Kit (Ready to Use) (orb552203)

Description	Human Relaxin/Insulin Like Family Peptide Receptor 3 (RXFP3) Sandwich ELISA KitReady to Use: Pre-dil...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NSD7
Sensitivity	0.057ng/mL
Entrez	51289
Expiration Date	Please enquire.