

Product Datasheet

Human SLAMF7 ELISA Kit (Ready to Use) (orb552385)

Description	Human Signaling Lymphocytic Activation Molecule Family, Member 7 (SLAMF7) Sandwich ELISA KitReady to...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NQ25
Sensitivity	0.053ng/mL
Entrez	57823
Expiration Date	Please enquire.