

Product Datasheet

Human SCCA1 ELISA Kit (Ready to Use) (orb552462)

Description	Human Squamous Cell Carcinoma Antigen 1 (SCCA1) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P29508
Sensitivity	34pg/mL
Entrez	6317
Expiration Date	Please enquire.