

Product Datasheet

Human STS ELISA Kit (Ready to Use) (orb552474)

Description	Human Steroid Sulfatase Isozyme S (STS) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	4.68-300ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P08842
Sensitivity	1.9ng/mL
Entrez	412
Expiration Date	Please enquire.