

Product Datasheet

Human STUB1 ELISA Kit (Ready to Use) (orb552477)

Description	Human STIP1 Homology And U-Box Containing Protein 1 (STUB1) Sandwich ELISA KitReady to Use: Pre-dilu...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.051ng/mL
Expiration Date	Please enquire.