

## Product Datasheet

# Human STIP1 ELISA Kit (Ready to Use) (orb552479)

<b>Description</b>	Human Stress Induced Phosphoprotein 1 (STIP1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ...
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P31948</a>
<b>Sensitivity</b>	6.5pg/mL
<b>Entrez</b>	<a href="#">10963</a>
<b>Expiration Date</b>	Please enquire.