

Product Datasheet

Human SPB ELISA Kit (Ready to Use) (orb552499)

Description	Human Surfactant Associated Protein B (SPB) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernatants and other biological fluids.
Uniprot ID	P07988
Sensitivity	0.54ng/mL
Entrez	6439
Expiration Date	Please enquire.