

Product Datasheet

Human SPD ELISA Kit (Ready to Use) (orb552501)

Description	Human Surfactant Associated Protein D (SPD) Sandwich ELISA Kit Ready to Use: Pre-diluted detection re...
Reactivity	Human
Range	6.25-400ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernatants and other biological fluids.
Uniprot ID	P35247
Sensitivity	2.16ng/mL
Entrez	6441
Expiration Date	Please enquire.