

Product Datasheet

Human SYP ELISA Kit (Ready to Use) (orb552507)

Description	Human Synaptophysin (SYP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P08247
Sensitivity	0.118ng/mL
Entrez	6855
Expiration Date	Please enquire.