

Product Datasheet

Human SYNPO ELISA Kit (Ready to Use) (orb552508)

| | |
|--------------------------|---|
| Description | Human Synaptopodin (SYNPO) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven... |
| Reactivity | Human |
| Range | 15.6-1000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids. |
| Sample Types | tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | Q8N3V7 |
| Sensitivity | 5.4pg/mL |
| Entrez | 11346 |
| Expiration Date | Please enquire. |