

Product Datasheet

Human TCIRG1 ELISA Kit (Ready to Use) (orb552526)

Description	Human T-Cell, Immune Regulator 1 (TCIRG1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum or other biological fluids.
Sample Types	serum and other biological fluids.
Uniprot ID	Q13488
Sensitivity	8.3pg/mL
Entrez	10312
Expiration Date	Please enquire.