

## Product Datasheet

### Human TAZ ELISA Kit (Ready to Use) (orb552528)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human Tafazzin (TAZ) Sandwich ELISA Kit<br>Ready to Use: Pre-diluted detection reagents for convenient a...         |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 0.156-10ng/mL   |
| <b>Form/Appearance</b>   | Ready to Use  |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.  |
| <b>Uniprot ID</b>        | <a href="#">Q16635</a>  |
| <b>Sensitivity</b>       | 0.054ng/mL  |
| <b>Entrez</b>            | <a href="#">6901</a>  |
| <b>Expiration Date</b>   | Please enquire.   |