

Product Datasheet

Human tPK1 ELISA Kit (Ready to Use) (orb552539)

Description	Human Tau-Protein Kinase 1 (tPK1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H3S4
Sensitivity	0.065ng/mL
Entrez	27010
Expiration Date	Please enquire.