

Product Datasheet

Human TNX ELISA Kit (Ready to Use) (orb552551)

Description	Human Tenascin X (TNX) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P22105
Sensitivity	0.53ng/mL
Entrez	7148
Expiration Date	Please enquire.