

Product Datasheet

Human TSPAN1 ELISA Kit (Ready to Use) (orb552558)

Description	Human Tetraspanin 1 (TSPAN1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conv...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O60635
Sensitivity	0.074ng/mL
Entrez	10103
Expiration Date	Please enquire.