

Product Datasheet

Human Tb4 ELISA Kit (Ready to Use) (orb552584)

Description	Human Thymosin Beta 4 (Tb4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P62328
Sensitivity	1.25ng/mL
Entrez	7114
Expiration Date	Please enquire.