

Product Datasheet

Human TRH ELISA Kit (Ready to Use) (orb552596)

Description	Human Thyrotropin Releasing Hormone (TRH) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	7.81-500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P20396
Sensitivity	3.25pg/mL
Entrez	7200
Expiration Date	Please enquire.