

Product Datasheet

Human TRAF3 ELISA Kit (Ready to Use) (orb552608)

Description	Human TNF Receptor Associated Factor 3 (TRAF3) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q13114
Sensitivity	29pg/mL
Entrez	7187
Expiration Date	Please enquire.