

Product Datasheet

Human TFF1 ELISA Kit (Ready to Use) (orb552669)

Description	Human Trefoil Factor 1 (TFF1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P04155
Sensitivity	54pg/mL
Entrez	7031
Expiration Date	Please enquire.