

Product Datasheet

Human TREM2 ELISA Kit (Ready to Use) (orb552675)

Description	Human Triggering Receptor Expressed On Myeloid Cells 2 (TREM2) Sandwich ELISA KitReady to Use: Pre-d...
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NZC2
Sensitivity	21.9pg/mL
Entrez	54209
Expiration Date	Please enquire.