

Product Datasheet

Human Tie1 ELISA Kit (Ready to Use) (orb552719)

| | |
|--------------------------|---|
| Description | Human Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) Sandwich ELISA ... |
| Reactivity | Human |
| Range | 1.56-100ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | P35590 |
| Sensitivity | 0.62ng/mL |
| Entrez | 7075 |
| Expiration Date | Please enquire. |