

Product Datasheet

Human VSIG2 ELISA Kit (Ready to Use) (orb552772)

Description	Human V-Set And Immunoglobulin Domain Containing Protein 2 (VSIG2) Sandwich ELISA KitReady to Use: P...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q96IQ7
Sensitivity	0.152ng/mL
Entrez	23584
Expiration Date	Please enquire.