

Product Datasheet

Human VS ELISA Kit (Ready to Use) (orb552791)

Description	Human Versican (VS) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient an...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q86W61
Sensitivity	27pg/mL
Entrez	1462
Expiration Date	Please enquire.