

Product Datasheet

Human vWFCP ELISA Kit (Ready to Use) (orb552817)

Description	Human Von Willebrand Factor Cleaving Protease (vWFCP) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q76LX8
Sensitivity	21.5pg/mL
Entrez	11093
Expiration Date	Please enquire.