

Product Datasheet

Human ZG16B ELISA Kit (Ready to Use) (orb552861)

Description	Human Zymogen Granule Protein 16 Homolog B (ZG16B) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q96DA0
Sensitivity	0.131ng/mL
Entrez	124220
Expiration Date	Please enquire.