

Product Datasheet

Mouse ART1 ELISA Kit (Ready to Use) (orb552895)

Description	Mouse ADP Ribosyltransferase 1 (ART1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
Reactivity	Mouse
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q60935
Sensitivity	0.32ng/mL
Entrez	11870
Expiration Date	Please enquire.