

Product Datasheet

Mouse a1M ELISA Kit (Ready to Use) (orb552923)

| | |
|--------------------------|---|
| Description | Mouse Alpha-1-Microglobulin/Bikunin Precursor (a1M) Sandwich ELISA KitReady to Use: Pre-diluted dete... |
| Reactivity | Mouse |
| Range | 1.56-100ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | Q07456 |
| Sensitivity | 0.52ng/mL |
| Entrez | 11699 |
| Expiration Date | Please enquire. |