

Product Datasheet

Mouse AIMP1 ELISA Kit (Ready to Use) (orb552934)

| | |
|--------------------------|--|
| Description | Mouse Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) Sandwich ELISA... |
| Reactivity | Mouse |
| Range | 31.2-2000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | Q3UZG4 |
| Sensitivity | 14.8pg/mL |
| Entrez | 13722 |
| Expiration Date | Please enquire. |