

Product Datasheet

Mouse Ab1-40 ELISA Kit (Ready to Use) (orb552938)

Description	Mouse Amyloid Beta Peptide 1-40 (Ab1-40) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	6.38pg/mL
Expiration Date	Please enquire.