

Product Datasheet

Mouse APPBP2 ELISA Kit (Ready to Use) (orb552942)

| | |
|--------------------------|---|
| Description | Mouse Amyloid Beta Precursor Protein Binding Protein 2 (APPBP2) Sandwich ELISA KitReady to Use: Pre-... |
| Reactivity | Mouse |
| Range | 0.312-20ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids. |
| Sample Types | tissue homogenates and other biological fluids. |
| Sensitivity | 0.115ng/mL |
| Expiration Date | Please enquire. |