

Product Datasheet

Mouse APP ELISA Kit (Ready to Use) (orb552945)

Description	Mouse Amyloid Precursor Protein (APP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
Reactivity	Mouse
Range	0.468-30ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P12023
Sensitivity	0.17ng/mL
Entrez	11820
Expiration Date	Please enquire.