

Product Datasheet

Mouse Angiotensinogen ELISA Kit (Ready to Use) (orb552948)

Description	Mouse Angiogenin (Angiotensinogen) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents fo...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P21570
Sensitivity	0.107ng/mL
Entrez	11727
Expiration Date	Please enquire.