

## Product Datasheet

# Mouse Angiotensinogen ELISA Kit (Ready to Use) (orb552957)

|                          |                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Angiostatin (Angiotensinogen)<br>Sandwich ELISA Kit<br>Ready to Use: Pre-diluted detection reagents f... |
| <b>Reactivity</b>        | Mouse                                                                                                          |
| <b>Range</b>             | 0.156-10ng/mL                                                                                                  |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                   |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                            |
| <b>Note</b>              | For research use only                                                                                          |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.                |
| <b>Sample Types</b>      | serum, plasma and other biological fluids.                                                                     |
| <b>Uniprot ID</b>        | <b>Q6TC10</b>                                                                                                  |
| <b>Sensitivity</b>       | 0.065ng/mL                                                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                                                                |