

Product Datasheet

Mouse ANXA5 ELISA Kit (Ready to Use) (orb552966)

Description	Mouse Annexin A5 (ANXA5) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P48036
Sensitivity	5.8pg/mL
Entrez	11747
Expiration Date	Please enquire.