

Product Datasheet

Mouse AMH ELISA Kit (Ready to Use) (orb552967)

Description	Mouse Anti-Mullerian Hormone (AMH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents fo...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q5EC55
Sensitivity	0.085ng/mL
Entrez	11705
Expiration Date	Please enquire.