

Product Datasheet

Mouse AP13 ELISA Kit (Ready to Use) (orb552970)

Description	Mouse Apelin 13 (AP13) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	3.9pg/mL
Expiration Date	Please enquire.