

Product Datasheet

Mouse CYP19A1 ELISA Kit (Ready to Use) (orb552999)

Description	Mouse Aromatase (CYP19A1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conveni...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.117ng/mL
Expiration Date	Please enquire.